

# PRECISION<sub>PAK</sub>

## Premeasured Concentrates

# GERMICIDAL CLEANER QT



### Technical Data

#### Where to Use

- Hospitals
- Child Care Centers
- Nursing Homes
- Restaurants
- Hotels/Motels
- Retail
- Health Care Facilities
- Transportation
- Schools
- Public Facilities
- Offices

#### Specifications

Item # .....84232-4  
 Size .....264 - 1 fl.oz. packs  
 Physical State .....Liquid  
 Color .....Light Red  
 Fragrance .....Fresh, Antiseptic  
 Specific Gravity .....1.00 +/- 0.007  
 Concentrate pH .....6  
 RTU pH .....10.2 +/- 0.6

HMIS HAZARD RATING - RTU solution	
HEALTH	2
FLAMMABILITY	0
REACTIVITY	0
PERSONAL PROTECTION	F

#### What the Pros Know

Germicidal Cleaner QT kills 99.9% of any disinfection organism listed on hard, non-porous surfaces. It is a one-step disinfectant that is effective against a broad spectrum of bacteria, is virucidal (including HIV-1, HCV & HBV) and inhibits the growth of mold and mildew and their odors when used as directed.

#### Description


- Disinfectant, cleaner, mildewstat, fungicide, virucide.
- Leaves fresh scent.
- Concentrated formula - each 1 oz. pack makes one quart of ready-to-use solution.

#### Directions

Wear PPE recommended on SDS, and read all directions and cautions prior to beginning your work.



Fill spray bottle with 32 oz. of water for spray clean.



Tear off one pack of SSS PrecisionPak Germicidal Cleaner QT. Carefully cut or tear corner of pack.



Mix in contents of one pack of Germicidal Cleaner QT into bottle. Rinse and discard pack.



For heavily soiled areas, a preliminary cleaning is required. Spray Germicidal Cleaner QT directly to soiled surfaces to thoroughly wet surfaces. Treated surfaces must remain wet for 10 minutes. Rinse with water and wipe off with a SSS MicroFiber Cloth (or comparable).

Just Add Water

Reduce Costs

Free Up Space

Leaves Fresh Scent

